**DB-2D and DB-3D**  
**Dri-Block®**  
Operator's Manual

Issue 23

Date of issue:

3rd May 2005.


Techne is a trade mark

© Techne 2002

**CONTENTS****DB-2D and DB-3D Operator's Manual**

	<i>page</i>
<b>SAFETY and INSTALATION</b>	
English	4
Français	6
Deutsch	8
Español	10
Italia	12
<b>THE DRI-BLOCK®</b>	14
<b>SPECIFICATION</b>	
DB-2D	15
DB-3D	16
Working Environment	17
Unpacking	17
<b>OPERATION</b>	18
Preparation	18
The overtemperature cutout	19
When you switch on	20
The front panel controls	20
Setting the Operating Temperature	22
After Use	23
<b>ADDITIONAL INFORMATION</b>	24
Insert Blocks	24
Operator maintenance	25
Calibration	26
Replacement parts	27

**NOTES**

## Introduction

Please read all the information in this booklet before using the unit.

## Warning

**HIGH TEMPERATURES ARE DANGEROUS:** they can cause serious burns to operators and ignite combustible material.

Techne have taken great care in the design of these units to protect operators from hazards, but Operators should pay attention to the following points:

- USE CARE AND WEAR PROTECTIVE GLOVES TO PROTECT HANDS;
- DO NOT put hot objects on or near combustible objects;
- DO NOT operate the unit close to inflammable liquids or gases;
- DO NOT place any liquid directly in your unit;
- At all times USE COMMON SENSE.

## Operator Safety

All Operators of Techne equipment must have available the relevant literature needed to ensure their safety. It is important that only suitably trained personnel operate this equipment, in accordance with the instructions contained in this manual and with general safety standards and procedures. If the equipment is used in a manner not specified by Techne the protection provided by the equipment to the Operator may be impaired.

All Techne units have been designed to conform to international safety requirements and are fitted with an overtemperature cutout. On some models, the cutout is adjustable and should be set to suit the application. On all other models the cutout is preset to protect the unit.

If a safety problem should be encountered, switch off at the mains socket and remove the plug from the supply.

## Installation

1. All Techne units are supplied with a power cable. This may be integral or plug-in.
2. Before connecting the mains supply, check the voltage against the rating plate. The rating plate is on the rear of the unit. Connect the mains cable to a suitable plug according to the table below.

**Note that the unit must be earthed to ensure proper electrical safety.**

<i>Connections</i>	<i>220V-240V</i>	<i>110V-120V</i>
Live	Brown	Black
Neutral	Blue	White
Earth	Green/yellow	Green

The fused plug supplied with the mains lead for use in the UK is fitted with the following value fuse to protect the cable: 5 AMP.

The fuse in the unit protects the unit and the operator

Note that units marked 230V on the rating plate work at 220V; units marked 120V work at 110V. In both cases, however, the heating rate will degrade by approximately 8%.

3. Plug the mains cable into the socket on the rear of the unit.
4. Place the unit on a suitable bench or flat workspace, or in a fume cupboard if required, ensuring that the air inlet vents on the underside are free from obstruction.
5. Note that the following symbols may be next to the indicator lamps on the front panel of the units and have the following meanings:

~ : the power indicator

⚡ : the heater indicator

🔥 : the overtemperature indicator

6. Symbols on or near the power switch of the unit have the following meanings:

I : mains switch On

O : mains switch Off

**After use**

When you have finished heating samples, remember that parts of the unit – the tubes, blocks and associated accessories – may be very hot. Take the precautions listed earlier.

**Guarantee**

The unit is guaranteed against any defect in material or workmanship for the period specified on the enclosed guarantee card. This period is from the date of purchase, and within this period all defective parts will be replaced free of charge provided that the defect is not the result of misuse, accident or negligence. Servicing under this guarantee should be obtained from the supplier.

Notwithstanding the description and specification(s) of the units contained in the Operator's Manual, Techne hereby reserves the right to make such changes as it sees fit to the units or to any component of the units.

This Manual has been prepared solely for the convenience of Techne customers and nothing in this Instruction Book shall be taken as a warranty, condition or representation concerning the description, merchantability, fitness for purpose or otherwise of the units or components.

**Operator maintenance**

NOTE: THAT THIS EQUIPMENT SHOULD ONLY BE DISMANTLED BY PROPERLY TRAINED PERSONNEL.

REMOVING THE SIDE, FRONT OR REAR PANELS EXPOSES POTENTIALLY LETHAL MAINS VOLTAGES.

THERE ARE NO OPERATOR MAINTAINABLE PARTS WITHIN THE EQUIPMENT.

In the unlikely event that you experience any problems with your unit which cannot easily be remedied, you should contact your supplier and return the unit if necessary. Please include any details of the fault observed and remember to return the unit in its original packing. Techne accept no responsibility for damage to units which are not properly packed for shipping: if in doubt, contact your supplier. See the Decontamination Certificate supplied with your unit.

**1. Cleaning**

Before cleaning your unit ALWAYS disconnect it from the power supply and allow it to cool below 50° C.

Your unit can be cleaned by wiping with a damp soapy cloth. Care should be exercised to prevent water from running inside the unit. Do not use abrasive cleaners.

**2. Overtemperature cutout**

In the event of no heater power, check the mains plug and lead. Repeated operation of the cutout indicates a serious fault: you may need to return the unit to your supplier for repair.

**3. Fuses**

Your unit is protected by one or two fuses. These should only be changed by suitably qualified personnel. If the fuses blow persistently, a serious fault is indicated and you may need to return the unit to your supplier for repair.

**Contact Information**

For technical, sales or servicing information, contact your local Techne dealer or,

Barloworld Scientific Ltd  
Beacon Road, Stone  
Staffordshire  
ST15 0SA, United Kingdom  
Telephone: +44(0)1785 812121  
Fax: +44(0)1785 813748  
e-mail: [equipment@barloworld-scientific.com](mailto:equipment@barloworld-scientific.com)

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## THE DRI-BLOCK®

Please read all the information in this manual before using the unit.

Techne's Dri-Block® heaters provide a safe, dry, constant temperature source in the laboratory. You can use the units for incubation, boiling, inactivation, wet ashing, sample concentration, enzyme analysis and many other clinical and industrial purposes.

The units cover the temperature range from 25°C to 200°C using machined aluminium alloy blocks as the heat transfer medium: temperature control circuits are built into the unit. Independent indicator lights indicate mains power, heater action and overtemperature cut-out. An ON/OFF switch is mounted on the rear panel, together with the mains connector and a double fuse holder.

The units are constructed in PPS plastic, which is strong yet lightweight and can resist sustained temperatures as high as 260°C. Even though the units heat up rapidly, highly efficient insulation ensures that cases remains cool enough to handle even at maximum operating temperatures.

Techne's Dri-Block® range has been designed to comply with all relevant RF (Radio Frequency) interference and electrical safety regulations.


The operating instructions in this Operator's Manual covers two of a range of Dri-Blocks®:

- DB-2D
- DB-3D

The features of individual models are outlined in the 'Specification' section of this manual.

Other units covered in other manuals are:

- DB-1M
- DB-2A
- DB-2P
- DB-3
- DB-3A


Typical rear connections

## SPECIFICATION

### The DB-2D Dri-Block®

This unit has a recessed chamber which can hold one or two interchangeable insert blocks which accept test tubes or other sample containers. The Temperature is set by 3 front panel buttons with digital display.

The DB-2D is supplied with an extraction tool for removing the insert blocks from the unit and two plastic half blocks for reducing heat losses when only one aluminium block is used.

#### Specification – DB-2D

Working temperature range	Ambient +5°C to 100°C
Settable temperature range	0.0°C to 105.0°C
Temperature stability	±0.05°C at 40°C
Temperature stability	±0.1°C at 100°C
Temperature setting	3 push buttons
Temperature display	4 digit LED
Set point resolution	0.1°C
Set point accuracy	±1°C
Max temperature variation between blocks	0.2°C at 40° C
Heat up time 30° C to 37° C	8 minutes
Heat up time 30° C to 56° C	9 minutes
Heat up time 30°C to 100°C	12 minutes
Number of blocks	2


#### Electrical supply

<i>Voltage</i>	<i>Cycles</i>	<i>Power</i>
230V	50Hz-60Hz	400W
120V	50Hz-60Hz	400W

## The DB-3D Dri-Block®

This unit has a recessed chamber which can hold up to three interchangeable insert blocks which accept test tubes or other sample containers. The Temperature is set by 3 front panel buttons with digital display.

The DB-3D is supplied with an extraction tool for removing the insert blocks from the unit.

Working temperature range	Ambient +5°C to 200°C
Settable temperature range	0.0°C to 200.0°C
Temperature stability	±0.1°C at 40°C
Temperature stability	±0.15°C at 100°C
Temperature setting	3 push buttons
Temperature display	4 digit LED
Set point resolution	0.1°C
Set point to accuracy	±1°C
Max temperature variation between blocks	0.2°C at 40° C
Heat up time 30° C to 37° C	11 minutes
Heat up time 30° C to 56° C	15 minutes
Heat up time 30°C to 200°C	25 minutes
Number of blocks	3

### Electrical supply

<i>Voltage</i>	<i>Cycles</i>	<i>Power</i>
230V	50Hz-60Hz	650W
110V-120V	50Hz-60Hz	650W

## Working Environment (all units)

The Dri-Block® units are designed to work safely under the following conditions:

Ambient temp. range 5°C to 40°C

Humidity Up to 95% relative humidity, non-condensing

**Note:** The control specifications are quoted for an ambient temperature of 20°C for units fitted with 2 or 3 insert blocks type F3506, or an F4467 in the case of the DB-1M. The specification will be closely held if the ambient temperature is in the range 10°C to 30°C. Outside this range the quoted figures may deteriorate but the unit will still work safely.

Radio frequency interference tested and passed to EN50081-1.

Immunity Tested and passed to EN50082-1

## Unpacking

When unpacking the unit, check that the following have been removed from the packing: The unit; a Guarantee card; a De-contamination Certificate.

Within the guarantee period, shown on the Guarantee Card, we undertake to supply replacements free of charge for parts which may on examination prove to be defective, provided that the defect is not the result of misuse, accident or negligence. Any instrument requiring service under this guarantee should be sent to the supplier through whom it was purchased, or, in the case of difficulty, it should be carefully packed in its original packing and consigned, carriage paid, to us. Techne takes no responsibility for returned goods damaged in transit.

**Returned goods will not be processed without a Returns Authorisation Number. Call + (44) (0)1785 812121 for a number.** On all correspondence, please quote the Serial Number in full and/or the Sales Order Number. **Please write the Returns Number on the outside of any packing.**

## OPERATION

### Preparation

1. Select the insert block(s) appropriate to the application. A list of available blocks is given on page 22 of this manual. Ensure that both the underside of the block(s) and the top of the hotplate are clean; efficient heat conduction between these two surfaces is essential.

**Top view showing three insert blocks in place**

2. Place the aluminium alloy insert blocks into the well of the unit and place the tubes containing the sample liquid in the blocks.
3. The heater design, temperature sensor and control circuit give good temperature control and uniformity, but make sure that there is a close fit of the tubes in the block to allow efficient heat transfer.
4. If less than the maximum number of insert blocks are fitted it is recommended that you fit plastic half blocks to the space either side of the aluminium block. The dummy blocks have the effect of reducing heat losses, thus improving temperature stability.
5. Plug the mains cable into the socket in the rear of the unit. Connect to the mains electricity supply with the plug provided or one wired correctly for your supply. Switch the power ON from the switch located next to the mains input lead. The mains indicator on the front panel will light (as well as the fuse neon on the mains input unit).

## Setting the Overtemperature Cut-out

### DB-2D

In the event of a fault occurring which could cause the heater to remain on, power is removed from the heater when the block temperature reaches approximately 110°C.

### DB-3D

An adjustable overtemperature cut-out is fitted to these units. As supplied it is set to remove power from the heater should the block temperature exceed approximately 210°C. If you need to prevent the block temperature exceeding a lower value, the overtemperature cut-out may be re-adjusted downwards. Note that the cut-out should be set to trip at least 5°C above the maximum desired operating temperature.

To set the overtemperature cut-out, turn the reset button (located on the back panel of the unit) fully clockwise then:

either:

- Heat the block to the desired cut-out temperature and turn the reset button anti-clockwise until the unit just cuts out.

or:


- Heat the block to the required maximum operating temperature and turn the reset button anti-clockwise until the unit just cuts out. Now turn the reset button an additional small amount clockwise.

When the cut-out operates, the unit stops working and the overtemperature cut-out indicator lights.

Remove the block from the unit by screwing the extractor tool into the thermometer hole and lifting the block vertically. Remember that high temperatures are dangerous; wear insulating gloves and take extreme care. Wait for the temperature to fall, perhaps by as much as 40°C, and press the reset button. The unit will now work as usual.

## When you Switch On

When you first switch on, the display will show the edition of the software which your unit uses. For example software issue "1" would be shown as follows:


It will display this for 1 second, then the actual temperature of the platten will be indicated.

## The Front Panel Controls


The front panel controls consist of three buttons for controlling the display, a four digit LED display and three indicators.

### The SET temperature Button


The **SET** temperature button displays the set temperature when pressed.

### The UP ARROW

When the **SET** temperature button is held down and the **UP ARROW** button is pressed, the set temperature is increased.

**The DOWN ARROW Button**

When the **SET** temperature button is held down and the **DOWN ARROW** button is pressed, the set temperature is decreased.

**Speed of Change of Set Temperature**


Each press of the **UP ARROW** or **DOWN ARROW** buttons will increase or decrease the set temperature by 0.1°C. If the buttons are held down the temperature change will accelerate to 5° per second

**Power Indicator**

The top indicator shows that there is power to the unit

**Heater Indicator**

The next indicator shows when the heater is heating the unit. When the temperature is being set, and the new set temperature is higher than the temperature already in the unit, the heater indicator will light as the unit tries to follow the set temperature. If the light is on continuously the heater is getting constant power. The only exception is described under Over-Temperature Indicator. As the temperature approaches the set temperature the heater indicator will flash. When set temperature is reached the indicator will stay on for shorter periods. If the Dri-Block is above the set temperature then the indicator will be off, as the heater is not getting any power.

**Over-Temperature Indicator**

If the unit should, for any reason, exceed the temperature set for the over-temperature cutout (see page 17), the over-temperature indicator will light. The heater will have been switched off and the unit will begin to cool even if the heater light is on (the light staying on or not depends on which circuit has sensed an over-temperature).


**Setting the Operating Temperature**

1. To display the set temperature on the digital display, press and hold the Set Temperature button.  
To adjust the set temperature, press the Set Temperature button and hold it while pressing the up or down buttons. When the Set Temperature button is released, the measured temperature is displayed (in degrees Celsius).
2. The heater (and heater indicator) comes on if the set temperature is higher than the current block temperature.
3. When the measured temperature approaches the set temperature, the heater indicator will begin to flash. As the measured temperature stabilises the indicator will stay on for shorter periods.
4. Due to variations in heat losses with different designs of insert block, the actual temperature may vary. The units are calibrated using model F3506 insert blocks. If you need to control the temperature to a greater accuracy than the instrument's set point, place a thermometer in the special hole in one of the blocks. For greater accuracy still, you can place the thermometer in the sample liquid in one of the test tubes. In either of these cases, it may be necessary to readjust the set temperature to achieve the precise temperature required. Allow the temperature of the unit to stabilise after each adjustment.

5. There will be a time lag between the heater platten and the insert block achieving the set temperature due to thermal contact between them. This is particularly true of the DB-3D where the display may reach temperature in 2 or 3 minutes while the insert block may take 12 minutes or more to reach temperature. Guide times are shown in the Specifications.

### After Use

1. When you have finished heating samples, remember that parts of the unit – the tubes, blocks and associated accessories – may be very hot. Take the precautions listed earlier. We recommend that the blocks should be allowed to cool to 70°C before being removed from the Dri-Block®. They will still have to be handled with care
2. If you need to remove an insert block while it is hot, screw the extractor tool into the thermometer hole and lift the block out vertically. Never leave the extractor tool in the block while it is being used in the Dri-Block®.

**Removing an insert block with the extractor tool**

**ADDITIONAL INFORMATION****Insert Blocks**

Insert blocks are made of aluminium alloy and must be ordered separately from the heater units. The following blocks are available from Techne:

<i>Part</i>	<i>Tube Size</i>	<i>Number of Holes</i>	<i>Size</i>
F3501	Plain block	None	95 x 76 x 51
F3502	6mm	30	95 x 76 x 51
F3503	10mm	20	95 x 76 x 51
F3504	12mm	20	95 x 76 x 51
F3505	13mm	20	95 x 76 x 51
F3506	15mm	12	95 x 76 x 51
F3507	16mm	12	95 x 76 x 51
F3508	19mm	8	95 x 76 x 51
F3509	25mm	6	95 x 76 x 51
F3510	10mm cuvettes	2 channels	95 x 76 x 51
F4460	Plain block	Thermometer hole only	95 x 76 x 51
F4461	9mm/7mm	10/20	95 x 76 x 51
F4462	24mm	6	95 x 76 x 51
F4463	26mm	6	95 x 76 x 51
F4464	1.5ml Eppendorf	20	95 x 76 x 51
F4465	0.5ml Eppendorf	30	95 x 76 x 51
F4466	Plastic half block	Plain	95 x 37 x 50
F4470	2 ml Eppendorf	20	95 x 76 x 51
F4471	0.2ml	72	95 x 76 x 51

## Operator Maintenance

NOTE THAT THIS EQUIPMENT SHOULD ONLY BE DISMANTLED BY PROPERLY TRAINED PERSONNEL. **REMOVING THE SIDE, FRONT OR REAR PANELS EXPOSES POTENTIALLY LETHAL MAINS VOLTAGES.** THERE ARE NO OPERATOR MAINTAINABLE PARTS WITHIN THE EQUIPMENT.

In the unlikely event that you experience any problems with your Dri-Block<sup>®</sup> which cannot easily be remedied, you should contact your supplier and return the unit if necessary. Please include any details of the fault observed and remember to return the unit in its original packing. Techne accept no responsibility for damage to units which are not properly packed for shipping: if in doubt, contact your supplier.

### 1. Cleaning

Before cleaning your unit ALWAYS disconnect from the power supply and allow to cool below 50° C.

Your Dri-Block<sup>®</sup> can be cleaned by wiping with a damp soapy cloth. Care should be exercised to prevent water from running inside the unit. Do not use abrasive cleaners.

Before using any cleaning or decontamination method except those recommended here, the responsible body or operator should check with Techne that the proposed method will not damage the equipment.

### 2. Overtemperature cut-out, DB-3D only

The overtemperature cut-out is a sensitive mechanical device and mechanical shock can cause it to trip.

- In the event of no heater power, check the mains plug and lead, then reset the cut-out control (located at the rear of the unit).
- Repeated operation of the cut-out indicates a serious fault: you may need to return the unit to your supplier for repair.

### 3. Fuses

Your unit is protected by two fuses. These should only be changed by suitably qualified personnel. If the fuses blow persistently, a serious fault is indicated and you may need to return the unit to your supplier for repair.

DB-2D	230V	2 x F2A;	120/100V	2 x F4A
DB-3D	230V	2 x F4A;	120/100V	2 x F6.3A

## Calibration

**Remember that if you change the calibration from that set at the factory you may change the calibration at all temperatures. You may get different calibration with different blocks.**

Make sure that you have good thermal contact between the blocks and the heater platen in the Dri-Block. In order to ensure that the calibration you are setting is correct, you will need to use an independent calibrated probe or thermometer.

Measure the actual temperature of the particular block you want to use in the Dri-Block using the calibrated probe or thermometer. If the calibration is not correct then you can follow this procedure. Make sure that the block and the unit are below 35°C.

- a Set the temperature display to 0.0°C. This will ensure that the block is at ambient temperature for the next procedure.
- b Hold down the Up and Down buttons and then press the Set button at the same time for 5 seconds. The display will change from the set temperature, 0.0, to the block temperature to indicate the start of self-calibration.
- c After a period of 17 minutes for the DB-2D or 25.5 minutes for the DB-3D the display will change from the block temperature (eg 40.0°C) to "EEEE". Measure the actual temperature of the block using the thermometer or probe. Press the 'SET' button and either the 'UP' or the 'DOWN' button to adjust the display to the same temperature as the measured value.
- d Press the 'UP' and 'DOWN' buttons together to confirm the value. The display will return to the block temperature and the heater will come full on.
- e After a period of 29.75 minutes for the DB-2D or 34 minutes for the DB-3D the display will change from the block temperature (eg 95.0°C for the DB-2D or 195°C for the DB-3D) to "EEEE".
- f Measure the actual temperature of the block using the thermometer. Press the 'SET' button and either the 'UP' or 'DOWN' button to adjust the display to the same temperature as the measured value.
- g Press the 'UP' and the 'DOWN' buttons together to confirm the value. The display will return to the block temperature and the unit will control with the new calibration parameters.

- h Set the temperature to the particular temperature you require. Measure the temperature when the block temperature has stabilised. Repeat the procedure as necessary; you may need to put in an off-set.

### Replacement Parts

The following parts may be obtained from your Techne dealer if replacements or alternatives are required:

<i>Part Number</i>	<i>Description</i>
6101308	Extraction tool
FCABLEUK	Mains cable with 13amp UK plug (5 amp fuse)
FCABLEEU	Mains cable with R/A Schuko plug
FCABLEUS	Mains cable with 3 pin US plug
FMW11	Techne HI-TEMP 96™ well plate, pack of 25
FMW12	HI-TEMP 96™ well cover, pack of 50
F4466	Plastic half block