CE350
Calibrator
Operator's Manual

Issue 4
Date of issue:
23rd July 2003.
Techne is the trade mark of
Techne and Techne Inc
© Techne 2003

CONTENTS

SAFETY and INSTALLATION	<i>page</i>
Declaration of Conformity	3
English	4
Français	6
Deutsch	8
Español	10
The CE350	12
Specification	13
Working environment	13
Warning	13
OPERATION	14
Preparation	14
Setting the operating temperature	14
After use	14
Changing from °C to °F	15
ADDITIONAL INFORMATION	16
Replacement parts	17
Alternative Blocks in Different Units	17

Issue 4 Date of issue: 23rd July 2003

Declaration of Conformity

Techne Unit CE350 has been designed to comply with the following European Standards:

EN 50081-1: 1992 Electromagnetic Compatibility - Generic emission standard.

EN 50082-1: 1992 Electromagnetic Compatibility - Generic immunity standard
(Performance criterion B).

EN 61010-1: 1993 Safety requirements for electrical equipment for measurement,
control and laboratory use.

Also:

IEC 1010-2-010: 1992 Particular requirements for laboratory equipment for the heating of
materials.

I have made all reasonable enquiries regarding the unit stated and its conformance to the
EU, Low Voltage and EMC Directives. To the best of my knowledge and belief this unit
conforms to these directives.

This Declaration is controlled under an ISO 9001:2000 system certificated by BSI Quality
Assurance, certificate number FM13585.

Signature	
Name	B C Coombes
Position	Quality Manager
Issue 3	23/07/2003

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Warning

HIGH TEMPERATURES ARE DANGEROUS: they can cause serious burns to operators and ignite combustible material.

Techne have taken great care in the design of these units to protect operators from hazards, but Operators should pay attention to the following points:

- USE CARE AND WEAR PROTECTIVE GLOVES TO PROTECT HANDS;
- DO NOT put hot objects on or near combustible objects;
- DO NOT operate the unit close to inflammable liquids or gases;
- DO NOT place any liquid directly in your unit;
- At all times USE COMMON SENSE.

Operator Safety

All Operators of Techne equipment must have available the relevant literature needed to ensure their safety. It is important that only suitably trained personnel operate this equipment, in accordance with the instructions contained in this manual and with general safety standards and procedures. If the equipment is used in a manner not specified by Techne the protection provided by the equipment to the Operator may be impaired.

All Techne units have been designed to conform to international safety requirements and are fitted with an over-temperature cutout. On some models, the cutout is adjustable and should be set to suit the application. On all other models the cutout is preset to protect the unit.

If a safety problem should be encountered, switch off at the mains socket and remove the plug from the supply.

Installation

1. All Techne units are supplied with a power cable.
2. Before connecting the mains supply, check the voltage against the rating plate. Connect the mains cable to a suitable plug according to the table below. Note that the unit must be earthed to ensure proper electrical safety.

<i>Connections 220V-240V</i>		<i>100V-120V</i>
Live	Brown	Black
Neutral	Blue	White
Earth	Green/yellow	Green

The fused plug supplied with the mains lead for use in the UK is fitted with the following value fuse to protect the cable: 230V UK 5 AMP

The fuse in the unit protects the unit and the operator

Note that units marked 230V on the rating plate work at 220V; units marked 120V work at 110V. In both cases, however, the heating rate will degrade by approximately 8%. The rating plate is on the rear of the unit.

3. Plug the mains cable into the socket on the rear of the unit.
4. Place the unit on a suitable bench or flat workspace, or in a fume cupboard if required, ensuring that the air inlet vents on the underside are free from obstruction.
5. Note that the following symbols may be next to the indicator lamps on the front panel of the units and have the following meanings:

 : the heater indicator

 : the over-temperature indicator

After use

When you have finished heating samples, remember that parts of the unit may be very hot. Take the precautions listed earlier.

Guarantee

The unit is guaranteed against any defect in material or workmanship for the period specified on the enclosed guarantee card. This period is from the date of purchase, and within this period all defective parts will be replaced free of charge provided that the defect is not the result of misuse, accident or negligence. Servicing under this guarantee should be obtained from the supplier.

Notwithstanding the description and specification(s) of the units contained in the Operator's Manual, Techne hereby reserves the right to make such changes as it sees fit to the units or to any component of the units.

This Manual has been prepared solely for the convenience of Techne customers and nothing in this Instruction Book shall be taken as a warranty, condition or representation concerning the description, merchantability, fitness for purpose or otherwise of the units or components.

Operator maintenance

NOTE THAT THIS EQUIPMENT SHOULD ONLY BE DISMANTLED BY PROPERLY TRAINED PERSONNEL. REMOVING THE SIDE, FRONT OR REAR PANELS EXPOSES POTENTIALLY LETHAL MAINS VOLTAGES. THERE ARE NO OPERATOR MAINTAINABLE PARTS WITHIN THE EQUIPMENT.

In the unlikely event that you experience any problems with your unit which cannot easily be remedied, you should contact your supplier and return the unit if necessary. Please include any details of the fault observed and remember to return the unit in its original packing. Techne accept no responsibility for damage to units which are not properly packed for shipping: if in doubt, contact your supplier. See the Decontamination Certificate supplied with your unit.

1. Cleaning

Before cleaning your unit ALWAYS disconnect it from the power supply and allow it to cool below 50° C.

Your unit can be cleaned by wiping with a damp soapy cloth. Care should be exercised to prevent water from running inside the unit. Do not use abrasive cleaners.

2. Over-temperature cutout

- In the event of no heater power, check the mains plug and lead.
- Repeated operation of the cutout indicates a serious fault: you may need to return the unit to your supplier for repair.

3. Fuses

Your unit is protected by one fuse. This should only be changed by suitably qualified personnel.

If the fuse blows persistently, a serious fault is indicated and you may need to return the unit to your supplier for repair.

Contact Information

For technical, sales or servicing information, contact your local Techne dealer or,
Techne, Duxford,

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THE CE 350

Please read all the information in this booklet before using the unit.

Techne's CE350 calibrator heater provides a safe, dry, constant temperature source. It is fast and economical and can be used either hand held or on a bench top. The weight of the unit is only 1.5 kilograms

The unit covers the temperature range from 20°C above ambient up to 350°C using machined aluminium alloy blocks as the heat transfer medium. The temperature control circuit is built into the unit.


Features include:

- Large block in 6 formats;
- Maximum temperature of 350°C;
- An independent over-temperature cutout;
- Temperature sensor burnout protection;
- Separate heater on/off switch for fast cool down without changing the set temperature;
- Indicators for over-temperature cutout and heater power.

Even though the unit heats up rapidly, highly efficient insulation ensures that the case remains cool enough to handle even at maximum operating temperatures.

Techne's calibrator has been designed to comply with all relevant Radio Frequency interference and electrical safety regulations.

Front Panel

Specification

Figures quoted are at the base of the well at the time of calibration.

This unit has a fixed block with six variations to hold four or five probes of different sizes.

See page 17 for the different blocks.

Temperature range	20 °C above ambient to 350 °C
Over-temperature limit	400 °C approximately
Temperature stability(10 min)	±0.15 °C @ 200°C
Display resolution	0.1 °C
Accuracy at 100°C to 300°C	±1°C
Heat up time 20° C to 300°C	5.2 minutes
Cool down 300°C to 100°C	10.8 minutes
Immersion Depth	101.4mm (4")
Fan Cooling	Automatic
Weight	1.5 Kg
Dimensions H x W x D	72 x 128 x 178 mm excluding the handle

Electrical supply

<i>Voltage</i>	<i>Cycles</i>	<i>Power</i>
230V	50/60Hz	400W
120V	50/60Hz	400W
100V	50/60Hz	400W

Working environment

The calibrator units are designed to work safely under the following conditions:

Ambient temp. range	5°C to 40°C
Humidity	Up to 95% relative humidity, non-condensing

Note: The control specifications are quoted for an ambient temperature of 20°C. The specification will be closely held if the ambient temperature is in the range 10°C to 30°C. Outside this range the quoted figures may deteriorate but the unit will still work safely.

Radio frequency interference tested and passed to EN50081-1.

Immunity Tested and passed to EN50082-1

Warning

Warning: HIGH TEMPERATURES ARE DANGEROUS
Avertissement: DANGER DE TEMPERATURES ELEVEES
Warnung: HOHE TEMPERATUREN SIND GEFÄHRLICH
Aviso: LAS TEMPERATURAS ELEVADAS SON PELI

OPERATION

Preparation

1. The heater design, temperature sensor and control circuit give good temperature control and uniformity, but make sure that there is a close fit of the probes in the block to allow efficient heat transfer.
3. Plug the power cable into the socket in the back of the unit. Connect the power cable to the electricity supply to the mains cable.

Setting the operating temperature

1. Ensure that the heater is switched off. The ON/OFF button should be in the out position not the latched position. The temperature can be set while the heater is on but the unit will try to follow the set temperature as it is changed.
2. To set the temperature required, press and hold the star button, *, while pressing the up arrow, ▲, or down arrow, ▼, button.
3. When you have the correct temperature displayed, release the star button, *, and the actual temperature of the block will be displayed.
4. When you are ready, press the heater ON/OFF button. The heater (and heater indicator) comes on if the set temperature is higher than the current block temperature. The temperature displayed will increase until it reaches the set temperature.
5. Allow the temperature to stabilise and calibrate your probe.
6. The controller is factory pre-set and only the set-temperature can be changed.

After use

1. When you have finished calibrating, remember that parts of the unit and the probe may be very hot. Take the precautions listed on page 4.
2. Switch off the heater by pressing and releasing the ON/OFF button. This will allow the fan to cool the block near to ambient.

Changing from °C to °F


If you wish to change the display units then you will have to follow this procedure:

Be careful not to change any of the other settings! The correct values are given below so that should you accidentally change one you can correct it.

1. Press both arrow buttons, ▲ ▼, for 3 sec. **tunE** will be displayed.
2. Press the down arrow button, ▼. **LEUL** will be displayed.
3. Press and hold the star button, *, and press the up arrow button, ▲, until **3** is displayed.
4. Press the up arrow button on its own, ▲, until **UEr** is displayed.
5. Press and hold both the arrow buttons, ▲ ▼, for 10 sec. **LoCK** will be displayed.
6. Press and hold the star button,*, and press the down arrow button, ▼, until **nonE** is displayed.
7. Press the down arrow button on its own, ▼, until **LEUL** is displayed.
8. Press and hold the star button,*, and press the down arrow button, ▼, until **2** is displayed.
9. Press and hold the up arrow button, ▲, until **unit** is displayed.
10. Press the star button,*, and the up or down arrow button, ▲ or ▼, to select °C or °F only.
11. When you have set the units you require:
Press the down arrow button, ▼, until **LEUL** is displayed.
12. Press and hold the star button, *, and press the up arrow button, ▲, until **3** is displayed.
13. Press the up arrow button, ▲, until **UEr** is displayed.
14. Press and hold both the arrow buttons, ▲ ▼, for 10 sec. **LoCK** will be displayed.
15. Press and hold the star button,*, and press the up arrow button, ▲, until **ALL** is displayed.
16. Press both the arrow buttons, ▲ ▼, for 3 sec to come out of the program mode.

Level 1	tune	off	Level 3	SP1.d	SSd
	bAnd	19		SP2.D	rly
	int.t	1.3		burn	uP.sc
	dEr.t	5		rEU.d	lr.2d
	dAC	3.0		rEU.L	ln.2n
	CYC.t	0.3		SPAn	0.0 #Adjust to suit
	ofSt	0		ZErø	0.0
	SP.LK	off		ChEK	off
	SEt.2	0		rEAd	UAr°
	bnd.2	1.0		dAtA	CtA
	CYC.2	on of		UEr	2
				rSet	none
Level 2	SP1.P	0	Level 4	LoCK	ALL
	hAnd	off		ProG	Auto
	PL.1	100		no.AI	off
	PL.2	100		diS.S	6
	SP2.A	none		dEr.S	0.5
	SP2.b	none			
	diSP	0.1°C			
	hi.SC	350.0			
	Lo.SC	0			
	inPt	Rtd (PRT)			
	unit	°C			

All units may be different, make a note of the value for your unit

ADDITIONAL INFORMATION

1. Before cleaning your unit ALWAYS disconnect from the power supply and allow to cool below 50° C.

Your calibrator can be cleaned by wiping with a damp soapy cloth. Care should be exercised to prevent water from running inside the unit. Do not use abrasive cleaners.

2. The controller is factory preset. See the Controller book for further details.

For safety reasons, most of the controller functions are locked. Please contact your dealer for further information

THE CONTROLLER MUST NOT BE SET TO CONTROL ABOVE 350°C OR DAMAGE TO THE UNIT MAY RESULT.

3. If the over-temperature cutout light comes on on the front panel, a fault has occurred and the cutout must reset before continuing. To reset the cutout press the red button on the base of the calibrator.

Repeated operation of the cutout indicates a serious fault: you may need to return the unit to your supplier for repair.

Replacement Parts

The following parts may be obtained from your Techne dealer if replacements or alternatives are required:


Techne Part Number Description

FCABLEUK	Mains cable with 13amp UK plug (5 amp fuse)
FCABLEEU	Mains cable with R/A Schuko plug
FCABLEUS	Mains cable with 3 pin US plug


Alternative Blocks in different units

Block Probe Sizes


A	4 x 6mm
B	1 each 10mm, 8mm, 6mm, 4.5mm, 3mm
C	2 each 6mm, 10mm
D	4 x 1/4"
E	1 each 3/8", 5/16", 1/4", 3/16", 1/8"
F	2 each 1/4", 3/8"

Blocks A and D

Blocks B and E

Blocks C and F