**DB-3 and DB-3A**  
**Dri-Block®**  
Operator's Manual

Issue 20

Date of issue:

23rd July 2003.

Techne is the trade mark of

Techne and Techne Inc

© Techne 2002

# CONTENTS

<b>DB-3 and DB-3A Operator's Manual</b>	<i>page</i>
<b>SAFETY and INSTALATION</b>	
CE Declaration of Conformity	3
English	4
Français	6
Deutsch	8
Español	10
Italia	12
<b>THE DRI-BLOCK®</b>	14
<b>SPECIFICATION</b>	
DB-3	15
DB-3A	16
Working Environment	17
<b>OPERATION</b>	18
Preparation	18
The overtemperature cutout	19
Setting the Operating Temperature	20
After Use	21
<b>ADDITIONAL INFORMATION</b>	22
Insert Blocks	22
Operator maintenance	23
Replacement parts	24

Techne  
Duxford  
Cambridge  
CB2 4PZ

### Declaration of Conformity

Techne Units DB-3 and DB-3A have been designed to comply with the following European Standards:

EN 50081-1:1992 Electromagnetic Compatibility; Generic emission standard.

EN 50082-1:1992 Electromagnetic Compatibility; Generic immunity standard (Performance criterion B).

EN 61010-1:1993 Safety requirements for electrical equipment for measurement, control and laboratory use.

EN 61010-2-010:1995 Particular requirements for laboratory equipment for the heating of materials.

I have made all reasonable enquiries regarding the unit stated and its conformance to the following EU directives.

Low Voltage directive, 73/23/EEC and amendment 93/68/EEC, and

EMC Directive 89/336/EEC and amendments 91/263/EEC 92/31/EEC and 93/68/EEC.

To the best of my knowledge and belief these units conform to these directives.

This Declaration is controlled under an ISO 9001:2000 system certificated by BSI Quality Assurance, certificate number FM13585.

Signature	<i>BCCoombes.</i>
Name	B C Coombes
Position	Quality Manager
Issue 5	23/07/2003

## Introduction

Please read all the information in this booklet before using the unit.

## Warning

**HIGH TEMPERATURES ARE DANGEROUS:** they can cause serious burns to operators and ignite combustible material.

Techne have taken great care in the design of these units to protect operators from hazards, but Operators should pay attention to the following points:

- USE CARE AND WEAR PROTECTIVE GLOVES TO PROTECT HANDS;
- DO NOT put hot objects on or near combustible objects;
- DO NOT operate the unit close to inflammable liquids or gases;
- DO NOT place any liquid directly in your unit;
- At all times USE COMMON SENSE.

## Operator Safety

All Operators of Techne equipment must have available the relevant literature needed to ensure their safety. It is important that only suitably trained personnel operate this equipment, in accordance with the instructions contained in this manual and with general safety standards and procedures. If the equipment is used in a manner not specified by Techne the protection provided by the equipment to the Operator may be impaired.

All Techne units have been designed to conform to international safety requirements and are fitted with an overtemperature cutout. On some models, the cutout is adjustable and should be set to suit the application. On all other models the cutout is preset to protect the unit.

If a safety problem should be encountered, switch off at the mains socket and remove the plug from the supply.

## Installation

1. All Techne units are supplied with a power cable. This may be integral or plug-in.
2. Before connecting the mains supply, check the voltage against the rating plate. The rating plate is on the rear of the unit. Connect the mains cable to a suitable plug according to the table below.

**Note that the unit must be earthed to ensure proper electrical safety.**

<i>Connections</i>	<i>220V-240V</i>	<i>110V-120V</i>
Live	Brown	Black
Neutral	Blue	White
Earth	Green/yellow	Green

The fused plug supplied with the mains lead for use in the UK is fitted with the following value fuse to protect the cable: 5 AMP.

The fuse in the unit protects the unit and the operator

Note that units marked 230V on the rating plate work at 220V; units marked 120V work at 110V. In both cases, however, the heating rate will degrade by approximately 8%.

3. Plug the mains cable into the socket on the rear of the unit.
4. Place the unit on a suitable bench or flat workspace, or in a fume cupboard if required, ensuring that the air inlet vents on the underside are free from obstruction.
5. Note that the following symbols may be next to the indicator lamps on the front panel of the units and have the following meanings:

 : the power indicator

 : the heater indicator

 : the overtemperature indicator

6. Symbols on or near the power switch of the unit have the following meanings:

I : mains switch On

O : mains switch Off

## After use

When you have finished heating samples, remember that parts of the unit – the tubes, blocks and associated accessories – may be very hot. Take the precautions listed earlier.

## Guarantee

The unit is guaranteed against any defect in material or workmanship for the period specified on the enclosed guarantee card. This period is from the date of purchase, and within this period all defective parts will be replaced free of charge provided that the defect is not the result of misuse, accident or negligence. Servicing under this guarantee should be obtained from the supplier.

Notwithstanding the description and specification(s) of the units contained in the Operator's Manual, Techne hereby reserves the right to make such changes as it sees fit to the units or to any component of the units.

This Manual has been prepared solely for the convenience of Techne customers and nothing in this Instruction Book shall be taken as a warranty, condition or representation concerning the description, merchantability, fitness for purpose or otherwise of the units or components.

## Operator maintenance

**NOTE: THAT THIS EQUIPMENT SHOULD ONLY BE DISMANTLED BY PROPERLY TRAINED PERSONNEL.**

**REMOVING THE SIDE, FRONT OR REAR PANELS EXPOSES POTENTIALLY LETHAL MAINS VOLTAGES.**

**THERE ARE NO OPERATOR MAINTAINABLE PARTS WITHIN THE EQUIPMENT.**

In the unlikely event that you experience any problems with your unit which cannot easily be remedied, you should contact your supplier and return the unit if necessary. Please include any details of the fault observed and remember to return the unit in its original packing. Techne accept no responsibility for damage to units which are not properly packed for shipping: if in doubt, contact your supplier. See the Decontamination Certificate supplied with your unit.

### 1. Cleaning

Before cleaning your unit ALWAYS disconnect it from the power supply and allow it to cool below 50° C.

Your unit can be cleaned by wiping with a damp soapy cloth. Care should be exercised to prevent water from running inside the unit. Do not use abrasive cleaners.

### 2. Overtemperature cutout

In the event of no heater power, check the mains plug and lead. Repeated operation of the cutout indicates a serious fault: you may need to return the unit to your supplier for repair.

### 3. Fuses

Your unit is protected by one or two fuses. These should only be changed by suitably qualified personnel. If the fuses blow persistently, a serious fault is indicated and you may need to return the unit to your supplier for repair.

## Contact Information

For technical, sales or servicing information, contact your local Techne dealer or,

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Duxford, CAMBRIDGE, CB2 4PZ, United Kingdom.

Telefono: +44 (0)1223 832401

Telefax: +44 (0)1223 836838

Assistenza: +44 (0)1223 836950

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e-mail: sales@technc.com

Web site: www.techne.com

## THE DRI-BLOCK®

Please read all the information in this manual before using the unit.

Techne's Dri-Block® heaters provide a safe, dry, constant temperature source in the laboratory. You can use the units for incubation, boiling, inactivation, wet ashing, sample concentration, enzyme analysis and many other clinical and industrial purposes.

The units cover the temperature range from 25°C to 200°C using machined aluminium alloy blocks as the heat transfer medium: temperature control circuits are built into the unit. Independent indicator lights indicate mains power, heater action and overtemperature cut-out. An ON/OFF switch is mounted on the rear panel, together with the mains connector and a double fuse holder.

The units are constructed in PPS plastic, which is strong yet lightweight and can resist sustained temperatures as high as 260°C. Even though the units heat up rapidly, highly efficient insulation ensures that cases remains cool enough to handle even at maximum operating temperatures.

Techne's Dri-Block® range has been designed to comply with all relevant RF (Radio Frequency) interference and electrical safety regulations.


The operating instructions in this Operator's Manual covers two of a range of Dri-Blocks®:

- DB-3
- DB-3A

The features of individual models are outlined in the 'Specification' section of this manual.

Other units covered in other manuals are:

- DB-1M
- DB-2A
- DB-2P
- DB-2D
- DB-3D


## SPECIFICATION

### The DB-3 Dri-Block®

This unit has a recessed chamber which can hold up to three interchangeable insert blocks which accept test tubes or other sample containers. Temperature is set by a front panel dial marked in degrees Celsius.

The DB-3 is supplied with an extraction tool for removing the insert blocks from the unit.

Temperature range

25°C to 100°C

Temperature stability

±0.1°C at 40°C

Temperature setting

Rotary knob

Temperature scale graduation 2°C

Set point accuracy ±2°C

Max temperature variation between blocks

0.2°C at 40°C

Heat up time 30°C to 37°C 8 minutes

Heat up time 30°C to 56°C 12 minutes

Heat up time 30°C to 100°C 18 minutes

Number of blocks 3


#### Electrical supply

<i>Voltage</i>	<i>Cycles</i>	<i>Power</i>
220V-240V	50Hz-60Hz	650W
110V-120V	50Hz-60Hz	650W
100V	50Hz-60Hz	650W

## The DB-3A Dri-Block®

This unit has a recessed chamber which can hold up to three interchangeable insert blocks which accept test tubes or other sample containers. Temperature is set by a front panel dial marked in degrees Celsius.

The DB-3A is supplied with an extraction tool for removing the insert blocks from the unit.

Temperature range	25°C to 200°C
Temperature stability	±0.1°C at 40°C
Temperature setting	Rotary knob
Temperature scale graduation	5°C
Set point accuracy	±2°C
Max temperature variation between blocks	0.2°C at 40° C
Heat up time 30° C to 37° C	8 minutes
Heat up time 30° C to 56° C	12 minutes
Heat up time 30°C to 200°C	30 minutes
Number of blocks	3

### Electrical supply

<i>Voltage</i>	<i>Cycles</i>	<i>Power</i>
220V-240V	50Hz-60Hz	650W
110V-120V	50Hz-60Hz	650W

## Working Environment (all units)

The Dri-Block® units are designed to work safely under the following conditions:

Ambient temp. range 5°C to 40°C

Humidity Up to 95% relative humidity, non-condensing

**Note:** The control specifications are quoted for an ambient temperature of 20°C for units fitted with 2 or 3 insert blocks type F3506, or an F4467 in the case of the DB-1M. The specification will be closely held if the ambient temperature is in the range 10°C to 30°C. Outside this range the quoted figures may deteriorate but the unit will still work safely.


Radio frequency interference tested and passed to EN50081-1.

Immunity Tested and passed to EN50082-1

## OPERATION

### Preparation

1. Select the insert block(s) appropriate to the application. A list of available blocks is given on page 22 of this manual. Ensure that both the underside of the block(s) and the top of the hotplate are clean; efficient heat conduction between these two surfaces is essential.

**Top view showing three insert blocks in place**

2. Place the aluminium alloy insert blocks into the well of the unit and place the tubes containing the sample liquid in the blocks.
3. The heater design, temperature sensor and control circuit give good temperature control and uniformity, but make sure that there is a close fit of the tubes in the block to allow efficient heat transfer.
4. If less than the maximum number of insert blocks are fitted it is recommended that you fit plastic half blocks to the space either side of the aluminium block. The dummy blocks have the effect of reducing heat losses, thus improving temperature stability.
5. Plug the mains cable into the socket in the rear of the unit. Connect to the mains electricity supply with the plug provided or one wired correctly for your supply. Switch the power ON from the switch located next to the mains input lead. The mains indicator on the front panel will light (as well as the fuse neon on the mains input unit).

## Setting the Overtemperature Cut-out

### DB-3

In the event of a fault occurring which could cause the heater to remain on, power is removed from the heater when the block temperature reaches approximately 110°C.

### DB-3A

An adjustable overtemperature cut-out is fitted to these units. As supplied it is set to remove power from the heater should the block temperature exceed approximately 210°C. If you need to prevent the block temperature exceeding a lower value, the overtemperature cut-out may be re-adjusted downwards. Note that the cut-out should be set to trip at least 5°C above the maximum desired operating temperature. To set the overtemperature cut-out, turn the reset button (located on the back panel of the unit) fully clockwise then:

either:

- Heat the block to the desired cut-out temperature and turn the reset button anti-clockwise until the unit just cuts out.

or:

- Heat the block to the required maximum operating temperature and turn the reset button anti-clockwise until the unit just cuts out. Now turn the reset button an additional small amount clockwise.

When the cut-out operates, the unit stops working and the overtemperature cut-out indicator lights.

Remove the block from the unit by screwing the extractor tool into the thermometer hole and lifting the block vertically. Remember that high temperatures are dangerous; wear insulating gloves and take extreme care. Wait for the temperature to fall, perhaps by as much as 40°C, and press the reset button. The unit will now work as usual.

## Setting the Operating Temperature


1. Turn the temperature control knob to the required temperature. Note that the actual temperature of the block critically depends on the type of block used and the quality of the thermal contact between the base of the block and the unit.
2. The heater (and heater indicator) comes on if the set temperature is higher than the current block temperature.
3. When the set temperature is approached, the heater indicator will begin to flash and when the temperature is stabilised the indicator will flash slowly.
4. Due to variations in heat losses with different designs of insert block, the actual temperature may vary. The units are calibrated using model F3506 insert blocks.

If you need to control the temperature to a greater accuracy than the instrument's set point, place a thermometer in the special hole in one of the blocks. For greater accuracy still, you can place the thermometer in the sample liquid in one of the test tubes. In either of these cases, it may be necessary to readjust the temperature control knob or preset adjustment to achieve the precise temperature required. Allow the temperature of the unit to stabilise after each adjustment.

5. There will be a time lag between the heater platten and the insert block achieving the set temperature due to thermal contact between them. The insert block may take up to 12 minutes to reach temperature, as shown in the Specifications on pages 15 and 16.

## After Use

1. When you have finished heating samples, remember that parts of the unit – the tubes, blocks and associated accessories – may be very hot. Take the precautions listed earlier. We recommend that the blocks should be allowed to cool to 70°C before being removed from the Dri-Block®. They will still have to be handled with care
2. If you need to remove an insert block while it is hot, screw the extractor tool into the thermometer hole and lift the block out vertically. Never leave the extractor tool in the block while it is being used in the Dri-Block®.

**Removing an insert block with the extractor tool**

## ADDITIONAL INFORMATION

### Insert Blocks

Insert blocks are made of aluminium alloy and must be ordered separately from the heater units. The following blocks are available from Techne:

<i>Part</i>	<i>Tube Size</i>	<i>Number of Holes</i>	<i>Size</i>
F3501	Plain block	None	95 x 76 x 51
F3502	6mm	30	95 x 76 x 51
F3503	10mm	20	95 x 76 x 51
F3504	12mm	20	95 x 76 x 51
F3505	13mm	20	95 x 76 x 51
F3506	15mm	12	95 x 76 x 51
F3507	16mm	12	95 x 76 x 51
F3508	19mm	8	95 x 76 x 51
F3509	25mm	6	95 x 76 x 51
F3510	10mm cuvettes	2 channels	95 x 76 x 51
F4460	Plain block	Thermometer hole only	95 x 76 x 51
F4461	9mm/7mm	10/20	95 x 76 x 51
F4462	24mm	6	95 x 76 x 51
F4463	26mm	6	95 x 76 x 51
F4464	1.5ml Eppendorf	20	95 x 76 x 51
F4465	0.5ml Eppendorf	30	95 x 76 x 51
F4466	Plastic half block	Plain	95 x 37 x 50
F4470	2 ml Eppendorf	20	95 x 76 x 51
F4471	0.2ml	72	95 x 76 x 51

## Operator Maintenance

NOTE THAT THIS EQUIPMENT SHOULD ONLY BE DISMANTLED BY PROPERLY TRAINED PERSONNEL. **REMOVING THE SIDE, FRONT OR REAR PANELS EXPOSES POTENTIALLY LETHAL MAINS VOLTAGES.** THERE ARE NO OPERATOR MAINTAINABLE PARTS WITHIN THE EQUIPMENT.

In the unlikely event that you experience any problems with your Dri-Block<sup>®</sup> which cannot easily be remedied, you should contact your supplier and return the unit if necessary. Please include any details of the fault observed and remember to return the unit in its original packing. Techne accept no responsibility for damage to units which are not properly packed for shipping: if in doubt, contact your supplier.

### 1. Cleaning

Before cleaning your unit ALWAYS disconnect from the power supply and allow to cool below 50° C.

Your Dri-Block<sup>®</sup> can be cleaned by wiping with a damp soapy cloth. Care should be exercised to prevent water from running inside the unit. Do not use abrasive cleaners.

Before using any cleaning or decontamination method except those recommended here, the responsible body or operator should check with Techne that the proposed method will not damage the equipment.

### 2. Overtemperature cut-out, DB-3A only

The overtemperature cut-out is a sensitive mechanical device and mechanical shock can cause it to trip.

- In the event of no heater power, check the mains plug and lead, then reset the cut-out control (located at the rear of the unit).
- Repeated operation of the cut-out indicates a serious fault: you may need to return the unit to your supplier for repair.

### 2. Fuses

Your unit is protected by two fuses. These should only be changed by suitably qualified personnel. If the fuses blow persistently, a serious fault is indicated and you may need to return the unit to your supplier for repair.

DB-3	230V	2 x F4A;	120/100V	2 x F6.3A
DB-3A	230V	2 x F4A;	120/100V	2 x F6.3A

## Replacement Parts

The following parts may be obtained from your Techne dealer if replacements or alternatives are required:

<i>Part Number</i>	<i>Description</i>
6101308	Extraction tool
FCABLEUK	Mains cable with 13amp UK plug (5 amp fuse)
FCABLEEU	Mains cable with R/A Schuko plug
FCABLEUS	Mains cable with 3 pin US plug
FMW11	Techne HI-TEMP 96™ well plate, pack of 25
FMW12	HI-TEMP 96™ well cover, pack of 50
F4466	Plastic half block