DC-3
Density Column
Clearing Device
Operator's Manual

Issue 11

Date of issue:

25th June 2003.

Techne is the trade mark of

Techne and Techne Inc

© Techne 2003

CONTENTS

DC-3 Operator's Manual	
	<i>page</i>
SAFETY and INSTALATION	
CE Declaration of Conformity	3
English	4
Français	6
Deutsch	8
Español	10
INTRODUCTION	12
Packing	12
Specification	12
Electrical Connection	13
OPERATION OF THE DC-3	14
SETUP OF THE DENSITY GRADIENT COLUMN	15
Liquid Systems	16
Liquid Preparation	17
Figure 1	18
Filling the Column	19
Calibrating the Column	19
INTRODUCING SAMPLES	20
Possible Sources of Error	21
Life of the Column	22
References	22
TECHNICAL INFORMATION	23
Operator maintenance	23
Replacement parts and Accessories	24

Techne
Duxford
Cambridge
CB2 4PZ

Declaration of Conformity

Techne Unit DC-3 has been designed to comply with the following European Standards:

EN 50081-1:1992 Electromagnetic Compatibility; Generic emission standard.

EN 50082-1:1992 Electromagnetic Compatibility; Generic immunity standard (Performance criterion B).

EN 61010-1:1993 Safety requirements for electrical equipment for measurement, control and laboratory use.

I have made all reasonable enquiries regarding the unit stated and its conformance to the following EU directives.

Low Voltage directive, 73/23/EEC and amendment 93/68/EEC, and

EMC Directive 89/336/EEC and amendments 91/263/EEC and 92/31/EEC .

To the best of my knowledge and belief these units conforms to these directives.

This Declaration is controlled under an ISO 9001:2000 system certificated by BSI Quality Assurance, certificate number FM13585.

Signature

Name

B C Coombes

Position

Quality Manager

Issue 5

23/07/2003

Introduction

Please read all the information in this booklet before using the unit.

Warning

HIGH TEMPERATURES ARE DANGEROUS: they can cause serious burns to operators and ignite combustible material.

Techne have taken great care in the design of these units to protect operators from hazards, but Operators should pay attention to the following points:

- USE CARE AND WEAR PROTECTIVE GLOVES TO PROTECT HANDS;
- DO NOT put hot objects on or near combustible objects;
- DO NOT operate the unit close to inflammable liquids or gases;
- DO NOT place any liquid directly in your unit;
- At all times USE COMMON SENSE.

Operator Safety

All Operators of Techne equipment must have available the relevant literature needed to ensure their safety.

It is important that only suitably trained personnel operate this equipment, in accordance with the instructions contained in this manual and with general safety standards and procedures. If the equipment is used in a manner not specified by Techne the protection provided by the equipment to the Operator may be impaired.

All Techne units have been designed to conform to international safety requirements and are fitted with an overtemperature cutout. On some models, the cutout is adjustable and should be set to suit the application. On all other models the cutout is preset to protect the unit.

If a safety problem should be encountered, switch off at the mains socket and remove the plug from the supply.

Installation

1. All Techne units are supplied with a power cable. This may be integral or plug-in.
2. Before connecting the mains supply, check the voltage against the rating plate. The rating plate is on the rear of the unit. Connect the mains cable to a suitable plug according to the table below.

Note that the unit must be earthed to ensure proper electrical safety.

<i>Connections</i>	<i>220V-240V</i>	<i>110V-120V</i>
Live	Brown	Black
Neutral	Blue	White
Earth	Green/yellow	Green

The fused plug supplied with the mains lead for use in the UK is fitted with the following value fuse to protect the cable: 10AMP.

The fuse in the unit protects the unit and the operator

Note that units marked 230V on the rating plate work at 220V; units marked 120V work at 110V. In both cases, however, the heating rate will degrade by approximately 8%.

3. Plug the mains cable into the socket on the rear of the unit.
4. Place the unit on a suitable bench or flat workspace, or in a fume cupboard if required, ensuring that the air inlet vents on the underside are free from obstruction.

After use

When you have finished heating samples, remember that parts of the unit – the tubes, blocks and associated accessories – may be very hot. Take the precautions listed earlier.

Guarantee

The unit is guaranteed against any defect in material or workmanship for the period specified on the enclosed guarantee card. This period is from the date of purchase, and within this period all defective parts will be replaced free of charge provided that the defect is not the result of misuse, accident or negligence. Servicing under this guarantee should be obtained from the supplier.

Notwithstanding the description and specification(s) of the units contained in the Operator's Manual, Techne hereby reserves the right to make such changes as it sees fit to the units or to any component of the units.

This Manual has been prepared solely for the convenience of Techne customers and nothing in this Instruction Book shall be taken as a warranty, condition or representation concerning the description, merchantability, fitness for purpose or otherwise of the units or components.

Operator maintenance

NOTE: THAT THIS EQUIPMENT SHOULD ONLY BE DISMANTLED BY PROPERLY TRAINED PERSONNEL.

REMOVING THE SIDE, FRONT OR REAR PANELS EXPOSES POTENTIALLY LETHAL MAINS VOLTAGES.

THERE ARE NO OPERATOR MAINTAINABLE PARTS WITHIN THE EQUIPMENT.

In the unlikely event that you experience any problems with your unit which cannot easily be remedied, you should contact your supplier and return the unit if necessary. Please include any details of the fault observed and remember to return the unit in its original packing. Techne accept no responsibility for damage to units which are not properly packed for shipping: if in doubt, contact your supplier. See the Decontamination Certificate supplied with your unit.

1.Cleaning

Before cleaning your unit ALWAYS disconnect it from the power supply and allow it to cool below 50° C.

Your unit can be cleaned by wiping with a damp soapy cloth. Care should be exercised to prevent water from running inside the unit. Do not use abrasive cleaners.

2.Fuses

Your unit is protected by one or two fuses. These should only be changed by suitably qualified personnel.

If the fuses blow persistently, a serious fault is indicated and you may need to return the unit to your supplier for repair.

Contact Information

For technical, sales or servicing information, contact your local Techne dealer or,
Techne, Duxford,
Cambridge, CB2 4PZ, United Kingdom.

Telephone:	+44(0)1223 832401
Fax:	+44(0)1223 836838
Service:	+44 (0)1223 836950 Out of office hours
e-mail:	sales@techne.com
Web site:	www.techne.com

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INTRODUCTION

Read the whole of this book before commencing work with the unit. Some of the DC-1 and DC-4 book is repeated here for your convenience. The clearing device consists of a conical plastic basket suspended by a nylon thread attached to a small reversible motor which is mounted on top of the plastic water jacket by two clamp screws. The basket is weighted to make it move axially within the graduated column tube and a fringe of polythene wipers around the edge of the basket ensures that all the specimens and floats within the column tube are collected as the basket is raised.

Packing

When you receive your unit, make sure you keep the original packing in case you ever need to return it for service or repair. The unit must be transported in the original packing to avoid damage. Techne accepts no responsibility for damage incurred unless the unit is correctly packed and transported in this way.

Specification

Housing	
Depth	114mm (4.5 in)
Width	94.4mm (3.6 in)
Height	114mm (4.5 in)
Basket	
Max diameter	50mm (1.9 in)
Length	54mm (2.1 in)
Thread	
Length when installed	914mm (36 in)
Material	Nylon
Speed	19.8mm/min (0.8 in/min)
Voltage	Cycles
220V/240V	50/60Hz
120V	50/60Hz
110V	50/60Hz

Electrical Connection

- 1 Before connecting the mains supply to the DC-3, check the voltage with that marked on the rating plate. Connect the mains cable to a suitable supply according to the table below. **Note that the unit must be earthed to ensure proper electrical safety.**

<i>CONNECTIONS</i>	<i>230V</i>	<i>110/120V</i>
Live	Brown	Black
Neutral	Blue	White
Earth	Green/yellow	Green

- 2 For connection of the Thermoregulator see the separate Operator's Manual.

OPERATION OF THE DC-3

- Set up the column.
- Wet the calibrated glass density floats and place them in the basket.
- Switch the motor to the 'DOWN' position.
- As each float reaches it's own density in the column it remains at that level while the basket continues to descend.
- The basket takes approximately 40 mins to travel the whole length of the column tube in order not to disturb the gradient.
- Switch the motor 'OFF' when the basket has reached the bottom of the tube.
- Make the column graph.
- The unit is ready to use.
- When it becomes necessary to clear the column tube of specimens, switch the motor to the 'UP' position and the basket will be drawn upwards collecting all the floats and specimens on the way.

SETUP OF THE DENSITY GRADIENT COLUMN

Jacket Filling

Tempette:

Replace the elbow connection from the base of the pump with the outlet pipe adaptor provided. Connect the plastic tube (inside the perspex vessel) to the adaptor.

- Mount the Tempette into the rectangular hole in the top plate of the water jacket. Lightly screw in the Tempette's clamp until it makes contact with the outside wall of the water jacket.
- Fill the vessel to within one inch of the top with demineralised water (through the column tube aperture) and then stand the whole unit on the floor (or on a bench if more convenient) near to where the column is to be utilised.

Note: A demineralised water with about 2% Ethylene Glycol type anti-freeze is a very useful alternative as it reduces deposits and algae growth, etc.

- A thermometer can be inserted into the vessel through the small hole provided in the top; set the Tempette (fully described in a separate User's Manual) dial to approx 23°C and allow the water to heat up. As the required temperature is approached, as shown on the thermometer, rotate the dial back until the neon indicator switches off at the exact temperature. Check the temperature is 23°C +0.1°C after the unit has been controlling for a short time.
- If the ambient temperature is higher than 19°C, then cold water should be circulated through the cooling coil in order that the Tempette can operate efficiently.
- Insert the column tube into the water vessel. To overcome the buoyancy, hold the column down by stretching the rubber band over the metal tags which are under the top cover fixing screws. The column is now ready for the filling operation.

Liquid Systems

The following list has been extracted from ASTM D.1505 - 60T.

System	Density Range <i>g/ml</i>
Methanol-Benzyl Alcohol	0.80 to 0.92
Isopropanol-Water	0.79 to 1.00
Isopropanol-Diethylene Glycol	0.79 to 1.11
Ethanol-Carbon Tetrachloride	0.79 to 1.59
Toluene-Carbon Tetrachloride	0.87 to 1.59
Water-Sodium Bromide	1.00 to 1.41
Water-Calcium Nitrate	1.00 to 1.60
Zinc Chloride-Ethanol-Water	0.80 to 1.70
Carbon Tetrachloride-1.3-Dibromopropane	1.60 to 1.99
1.3 Dibromopropane-Ethylene Bromide	1.99 to 2.18
Ethylene Bromide-Bromoform	2.18 to 2.89
Carbon Tetrachloride-Bromoform	1.60 to 2.89
Tetrachloroethylene-Bromoform	1.55 to 2.70

Liquid Preparation

- When selecting the two liquids, the following are important requirements:
 - 1 Both liquids should be readily available and they must obviously be miscible.
 - 2 Preferably colourless, low viscosity liquids.
 - 3 The liquids must obviously span the required density gradient, and should easily wet the samples.
 - 4 The liquids must not have an adverse effect on the samples to be tested.
 - 5 The two liquids can themselves be prepared solutions prior to mixing in the filling device.
 - Assemble the flasks as shown in the diagram.
 - Close the outlet taps on the two conical flasks.
 - With the above considerations in mind, prepare two liquids;
 - 'a' having a density 20% (of the whole range) below the lowest density required and
 - 'b' having a density 30% (of the whole range) above the highest density required (ref BS3715).
 - Place them in the separate conical flasks.
 - Where the densities of the two liquids differ by 5% or more, it is necessary to put a greater volume of the liquid of the lower density in its flask.
 - The heights of the surfaces of the two liquids above the filling device platform should be in inverse proportion to their density. This will ensure that when the connecting tap between the two is opened, an approximate equilibrium condition will exist.
 - It is important that the total volume of the two liquids should be two litres.
- Note:** The wider the column range required, the more difficult it becomes to obtain linearity for ranges not exceeding approximately 0.3 g/ml but with ranges exceeding this value it may be necessary to either:
- 1 Use an auxiliary stirrer through the top of the flask with the magnetic stirrer. In this instance care must be taken not to introduce air bubbles into the liquid. or
 - 2 Considerably increase the total filling time; up to, say, four hours.
- If necessary, the liquids may be degassed by evacuation at a pump pressure of about 60cm. A water filter pump is suitable.


Figure 1

Filling the Column

- Insert the capillary filling tube (item 2) down into the column tube until its lower end rests on the bottom.
- Link up the rest of the apparatus as shown in fig 1 using a small amount of silicon grease on all the joints; the stirrer and flasks can stand on a bench alongside the column.
- Connect the magnetic stirrer to the correct electrical supply and switch on.
- Fully open the connecting tap between the two vessels. Allow time for hydrostatic equilibrium to be established between the two vessels and then open the outlet tap to the capillary tube.
- The column will slowly fill, and should take about 2 hours.
- Regulation of filling rate is possible by adjustment of the outlet tap but this **must not** be done once a filling operation has started.
- During a filling operation the density of liquid 'a' will increase continuously as the liquid of higher density in flask 'b' flows into it. This will introduce a solution of continually graded density into the column.
- Shut off the outlet tap when the column tube is sufficiently filled above the scale top mark.
- Switch off the stirrer.
- Very slowly and carefully remove the capillary filling tube.
- Retain some of the denser liquid for subsequent specimen wetting.

Calibrating the Column

- Position the apparatus where it is to be used.
- Select eight calibrated glass density floats to cover the density range of the column, and ensure they are very clean.
- Starting with the highest value float, wet each one in some of the left over solution and place them into the column. The DC-3 may be used for this. Before actually releasing the floats ensure that they have no visible air bubbles attached.
- The floats will then assume their own density level in the column.
- By measuring the height of each float in the column and knowing the density of the float a calibration graph can be drawn.

If the readings obtained from several specimens of the same sample are widely different, then more specimens should be tested to see if the spread of results is real.

INTRODUCING SAMPLES

Since a large number of samples may eventually be placed into the column, some form of identification is necessary. This is best achieved by identifying the shape of the specimens in simple geometric forms. It is important that their centre of volume is easily known, and by combinations of geometric shapes a large number of different specimens is possible.

Always use a sharp cutting tool to avoid straining the material.

- After the specimens have reached equilibrium their positions (centres of volume) can be determined by means of the antiparallax scale on the column tube.
- The density of the specimen can then be read directly off the graph.
- When attempting to assess density to ± 0.0002 g/ml the specimens should always be double checked, even by immersion in a second column.
- The graph should be replotted each time the column is used.
- If practical it is advisable to allow at least 12 hours before taking readings of specimen density.

- 1 The specimens may have **voids** or have **air bubbles** attached.
In such cases the highest density value should be accepted especially if it is the most frequent one.
- 2 **Foreign matter** embedded in the specimens.
In this case, cut the specimen in to two ore more pieces and determine the density of these by replacing in a column.
- 3 The column becoming so **choked with specimens** that interference between them may prevent precise readings.
It is possible to clear the column of all specimens and floats by means of a mesh basket, DC-3, which is drawn very slowly up through the column. However, it is often economic to simply scrap the column and make a fresh one.

Life of the Column

- The life of a column varies greatly but normally lasts around fifteen weeks if undisturbed.
- However, sometimes a column might have to be discarded as soon as it is made because the graph obtained is unacceptable. This occurs particularly if there has been unevenness in the flow during filling. For example, if the capillary becomes blocked or if there is an air bubble in the capillary.
- It is impossible to give a quantitative guide to the expected life of a column because it depends on many factors, some of which may be impossible to control.
- In general it can be said that any disturbance that produces convection currents will hasten the deterioration compared with that calculated from the diffusion equation.
- Thus the service life will depend on the number of samples tested, the way the samples are inserted, the amount of differential evaporation, and temperature gradients.

Note: A column clearing device is available, but has been omitted as standard equipment because it is only necessary where either the rate of sample testing is high or if the column solutions are expensive.

References

'Measuring the density of polyolefins with an improved Gradient column' N. Payne and C.E. Stephenson. Materials Research & standards 4 3-7. (January 1964)

BS 2782:Part 6:Methods 620A to 620D:1991

ASTM D. 1505 - 60T.

TECHNICAL INFORMATION

Brief fault finding notes are given in this section.

NOTE THAT ELECTRICAL EQUIPMENT SHOULD ONLY BE DISMANTLED BY PROPERLY TRAINED PERSONNEL. REMOVING THE CASE EXPOSES POTENTIALLY LETHAL MAINS VOLTAGES. THERE ARE NO USER MAINTAINABLE PARTS WITHIN THIS EQUIPMENT.

User Maintenance

General

In the unlikely event that you experience any problems with your Density Column or accessories which cannot be easily remedied, you should contact your supplier and return the unit if necessary. Please include details of the fault observed and remember to return the unit in its original packing. Techne accept no responsibility for damage to units which are not properly packed for shipping: if in doubt, contact your supplier.

Fuses

If the power lamp, on the front panel, fails to light when power is supplied to the unit the fuses may have blown. Check that there is no external cause such as a faulty plug.

These should only be changed by a suitably qualified person. If the fuses blow persistently, a serious fault is indicated and you may need to return the unit to your supplier for repair.

Cleaning your Density Column

Before cleaning your unit ALWAYS disconnect from the power supply and allow to cool.

You can clean the case of the Filling Device with a cloth dampened in water or ethanol. Only use a damp cloth for the Density Column. No part of the unit should be immersed in solvents. Do not use acetone or abrasive cleaners.

Before using any cleaning or decontamination method except those recommended in this manual, the responsible body should check with Techne that the proposed method will not damage the equipment.

Replacement Parts and Accessories

DC-1 and DC-4

When ordering spare parts check fig 1, page 18, to obtain the correct item reference and always quote the instrument serial number.

<i>Item N°</i>	<i>Description</i>	<i>Qty</i>	<i>Notes</i>
1	F2019 CONICAL FLASK	2	
2	F2020 CAPILLARY TUBE	1	
3	F2021 STOPCOCK	1	
4	F2022 STOPCOCK	1	
5	F2023 JOINT CLIP	5	
6	F2024 BALL BEARING	1	
7	F2025 MAGNETIC STIRRER BAR	1	
8	5000821 9mm x 6mm PORTEX TUBE	1	
9	F1947 COOLING COIL	1	
10	6005730 REDUCING NOZZLE	1	
12	6009689 OUTLET PIPE ADAPTOR	1	
14	6002848 CABLE CLIP	2	
15	6100255 WATER JACKET	1	DC-1
16	F2015 GLASS COLUMN ASSEMBLY	1	DC-1
15	6101029 WATER JACKET	1	DC-4
16	F2015 GLASS COLUMN ASSEMBLY	2	DC-4
17	6002014 RUBBER CAP WITH HOLE	2	
17	6002027 RUBBER CAP	2	
18	5000823 5/16 x 3/16 RUBBER TUBING N.5	1	
34	FJP8D TE-8J THERMOREGULATOR	1	240V
34	FJP8P TE-8J THERMOREGULATOR	1	120V
34	FJP8Y TE-8J THERMOREGULATOR	1	100V

DC-3

<i>Number</i>	<i>Description</i>	<i>Qty</i>
1	6004806 BOBBIN	1
2	6004919 BASKET ASSEMBLY	1
3	6005924 LINE GUIDE	1
4	6001259 SHOULDER SCREW	2