**DB-1M, DB-2P
and DB-2A Dri-Block®
Operator's Manual**

Issue 21

Date of issue:

23rd July 2003.

Techne is the trade mark of

Techne and Techne Inc

© Techne 2002


Barloworld Scientific

CONTENTS

DB-1M, DB-2P and DB-2A Operator's Manual

	<i>page</i>
SAFETY and INSTALATION	
CE Declaration of Conformity	3
English	4
Français	6
Deutsch	8
Español	10
Italia	12
THE DRI-BLOCK®	14
Packing	15
Working Environment	15
SPECIFICATIONS	
DB-1M	16
DB-2P	17
DB-2A	18
OPERATION	19
Preparation	19
The overtemperature cutout	20
Setting the Pre-set Temperature	21
Setting the Operating Temperature	23
After Use	24
ADDITIONAL INFORMATION	25
Insert Blocks	25
Operator maintenance	26
Replacement parts	27

Techne
Duxford
Cambridge
CB2 4PZ

Declaration of Conformity

Techne Units DB-1M, DB-2P and DB-2A have been designed to comply with the following European Standards:

EN 50081-1:1992 Electromagnetic Compatibility; Generic emission standard.

EN 50082-1:1992 Electromagnetic Compatibility; Generic immunity standard (Performance criterion B).

EN 61010-1:1993 Safety requirements for electrical equipment for measurement, control and laboratory use.

EN 61010-2-010:1995 Particular requirements for laboratory equipment for the heating of materials.

I have made all reasonable enquiries regarding the unit stated and its conformance to the following EU directives.

Low Voltage directive, 73/23/EEC and amendment 93/68/EEC, and

EMC Directive 89/336/EEC and amendment 92/31/EEC.

To the best of my knowledge and belief these units conform to these directives.

This Declaration is controlled under an ISO 9001:2000 system certificated by BSI Quality Assurance, certificate number FM13585.

Signature

Name

B C Coombes

Position

Quality Manager

Issue 4

23/07/2003

Introduction

Please read all the information in this booklet before using the unit.

Warning

HIGH TEMPERATURES ARE DANGEROUS: they can cause serious burns to operators and ignite combustible material.

Techne have taken great care in the design of these units to protect operators from hazards, but Operators should pay attention to the following points:

- USE CARE AND WEAR PROTECTIVE GLOVES TO PROTECT HANDS;
- DO NOT put hot objects on or near combustible objects;
- DO NOT operate the unit close to inflammable liquids or gases;
- DO NOT place any liquid directly in your unit;
- At all times USE COMMON SENSE.

Operator Safety

All Operators of Techne equipment must have available the relevant literature needed to ensure their safety. It is important that only suitably trained personnel operate this equipment, in accordance with the instructions contained in this manual and with general safety standards and procedures. If the equipment is used in a manner not specified by Techne the protection provided by the equipment to the Operator may be impaired.

All Techne units have been designed to conform to international safety requirements and are fitted with an overtemperature cutout. On some models, the cutout is adjustable and should be set to suit the application. On all other models the cutout is preset to protect the unit.

If a safety problem should be encountered, switch off at the mains socket and remove the plug from the supply.

Installation

1. All Techne units are supplied with a power cable. This may be integral or plug-in.
2. Before connecting the mains supply, check the voltage against the rating plate. The rating plate is on the rear of the unit. Connect the mains cable to a suitable plug according to the table below.

Note that the unit must be earthed to ensure proper electrical safety.


<i>Connections</i>	<i>220V-240V</i>	<i>110V-120V</i>
Live	Brown	Black
Neutral	Blue	White
Earth	Green/yellow	Green

The fused plug supplied with the mains lead for use in the UK is fitted with the following value fuse to protect the cable: 5 AMP.

The fuse in the unit protects the unit and the operator

Note that units marked 230V on the rating plate work at 220V; units marked 120V work at 110V. In both cases, however, the heating rate will degrade by approximately 8%.

3. Plug the mains cable into the socket on the rear of the unit.
4. Place the unit on a suitable bench or flat workspace, or in a fume cupboard if required, ensuring that the air inlet vents on the underside are free from obstruction.
5. Note that the following symbols may be next to the indicator lamps on the front panel of the units and have the following meanings:

	:	the power indicator
	:	the heater indicator
	:	the overtemperature indicator

6. Symbols on or near the power switch of the unit have the following meanings:

I	:	mains switch On
O	:	mains switch Off

After use

When you have finished heating samples, remember that parts of the unit – the tubes, blocks and associated accessories – may be very hot. Take the precautions listed earlier.

Guarantee

The unit is guaranteed against any defect in material or workmanship for the period specified on the enclosed guarantee card. This period is from the date of purchase, and within this period all defective parts will be replaced free of charge provided that the defect is not the result of misuse, accident or negligence. Servicing under this guarantee should be obtained from the supplier.

Notwithstanding the description and specification(s) of the units contained in the Operator's Manual, Techne hereby reserves the right to make such changes as it sees fit to the units or to any component of the units.

This Manual has been prepared solely for the convenience of Techne customers and nothing in this Instruction Book shall be taken as a warranty, condition or representation concerning the description, merchantability, fitness for purpose or otherwise of the units or components.

Operator maintenance

NOTE: THAT THIS EQUIPMENT SHOULD ONLY BE DISMANTLED BY PROPERLY TRAINED PERSONNEL.

REMOVING THE SIDE, FRONT OR REAR PANELS EXPOSES POTENTIALLY LETHAL MAINS VOLTAGES.

THERE ARE NO OPERATOR MAINTAINABLE PARTS WITHIN THE EQUIPMENT.

In the unlikely event that you experience any problems with your unit which cannot easily be remedied, you should contact your supplier and return the unit if necessary. Please include any details of the fault observed and remember to return the unit in its original packing. Techne accept no responsibility for damage to units which are not properly packed for shipping; if in doubt, contact your supplier. See the Decontamination Certificate supplied with your unit.

1. Cleaning

Before cleaning your unit ALWAYS disconnect it from the power supply and allow it to cool below 50° C.

Your unit can be cleaned by wiping with a damp soapy cloth. Care should be exercised to prevent water from running inside the unit. Do not use abrasive cleaners.

2. Overtemperature cutout

In the event of no heater power, check the mains plug and lead. Repeated operation of the cutout indicates a serious fault: you may need to return the unit to your supplier for repair.

3. Fuses

Your unit is protected by one or two fuses. These should only be changed by suitably qualified personnel. If the fuses blow persistently, a serious fault is indicated and you may need to return the unit to your supplier for repair.

Contact Information

For technical, sales or servicing information, contact your local Techne dealer or,

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Duxford, CAMBRIDGE, CB2 4PZ, United Kingdom.

Telefono: +44 (0)1223 832401

Telefax: +44 (0)1223 836838

Assistenza: +44 (0)1223 836950

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e-mail: sales@technc.com

Web site: www.techne.com

THE DRI-BLOCK®

Please read all the information in this manual before using the unit.

Techne's Dri-Block® heaters provide a safe, dry, constant temperature source in the laboratory. You can use the units for incubation, boiling, inactivation, wet ashing, sample concentration, enzyme analysis and many other clinical and industrial purposes.

The units cover the temperature range from 25°C to 200°C using machined aluminium alloy blocks as the heat transfer medium: temperature control circuits are built into the unit. Independent indicator lights indicate mains power, heater action and overtemperature cut-out. An ON/OFF switch is mounted on the rear panel, together with the mains connector and a double fuse holder.

The units are constructed in PPS plastic, which is strong yet lightweight and can resist sustained temperatures as high as 260°C. Even though the units heat up rapidly, highly efficient insulation ensures that cases remain cool enough to handle even at maximum operating temperatures.

Techne's Dri-Block® range has been designed to comply with all relevant RF (Radio Frequency) interference and electrical safety regulations.


The operating instructions in this Operator's Manual covers three of a range of Dri-Blocks®:

- DB-1M
- DB-2A
- DB-2P

The features of individual models are outlined in the 'Specification' section of this manual.

Other units covered in other manuals are:

- DB-3
- DB-3A
- DB-2D
- DB-3D

Unpacking

When unpacking the unit, check that the following have been removed from the packing: The unit; a Guarantee card; a De-contamination Certificate.

Within the guarantee period, shown on the Guarantee Card, we undertake to supply replacements free of charge for parts which may on examination prove to be defective, provided that the defect is not the result of misuse, accident or negligence. Any instrument requiring service under this guarantee should be sent to the supplier through whom it was purchased, or, in the case of difficulty, it should be carefully packed in its original packing and consigned, carriage paid, to us. Techne takes no responsibility for returned goods damaged in transit.

Returned goods will not be processed without a Returns Authorisation Number. Call + (44) (0)1223 832401 for a number. On all correspondence, please quote the Serial Number in full and/or the Sales Order Number. Please write the Returns Number on the outside of any packing.

Working Environment (all units)

The Dri-Block® units are designed to work safely under the following conditions:

Ambient temp. range 5°C to 40°C

Humidity Up to 95% relative humidity, non-condensing

Note: The control specifications are quoted for an ambient temperature of 20°C for units fitted with 2 or 3 insert blocks type F3506, or an F4467 in the case of the DB-1M. The specification will be closely held if the ambient temperature is in the range 10°C to 30°C. Outside this range the quoted figures may deteriorate but the unit will still work safely.

Radio frequency interference tested and passed to EN50081-1.


Immunity Tested and passed to EN50082-1

SPECIFICATION

The DB-1M Dri-Block®

This unit has a recessed chamber to hold interchangeable blocks. The blocks are machined to accept various microplates, one of which is the Techne HI-TEMP 96™ microplate. Temperature control is provided by three adjustable presets on the front panel. The DB-1M has a cover to hold the plates in optimum thermal contact with the block.

The DB-1M is supplied with a screw driver for adjusting the pre-sets and a plastic lid which clips onto the hinge on the unit.

Temperature range	25°C to 105°C
Temperature stability	±0.1°C at 40°C
Temperature setting	3 adjustable presets
Max. temperature variation	0.2°C at 40° C
Heat up time 30° C to 37° C	15 minutes
Heat up time 30° C to 56° C	27 minutes
Heat up time 30° C to 100° C	30 minutes
Number of blocks	1 Multiwell


Electrical supply

<i>Voltage</i>	<i>Cycles</i>	<i>Power</i>
220V-240V	50Hz-60Hz	400W
110V-120V	50Hz-60Hz	400W

The DB-2P Dri-Block®

This unit has a recessed chamber to hold one or two interchangeable insert blocks which accept test tubes or other sample containers. Temperature is set by three adjustable presets on the front panel.

The DB-2P is supplied with a screw driver for adjusting the pre-sets, an extraction tool for removing the insert blocks from the unit and two plastic half blocks for reducing heat losses when only one aluminium block is used.

Temperature range	25°C to 105°C
Temperature stability	±0.1°C at 40°C
Temperature setting	3 adjustable presets
Max temperature variation between blocks	0.2°C at 40° C
Heat up time 30° C to 37° C	9 minutes
Heat up time 30° C to 56° C	10 minutes
Heat up time 30°C to 100°C	15 minutes
Number of blocks	2


Electrical supply

<i>Voltage</i>	<i>Cycles</i>	<i>Power</i>
220V-240V	50Hz-60Hz	400W
110V-120V	50Hz-60Hz	400W

The DB-2A Dri-Block®

This unit has a recessed chamber which can hold one or two interchangeable insert blocks which accept test tubes or other sample containers. Temperature is set by a front panel dial marked in degrees Celsius.

The DB-2A is supplied with an extraction tool for removing the insert blocks from the unit and two plastic half blocks for reducing heat losses when only one aluminium block is used.

Temperature range	25°C to 100°C
Temperature stability	±0.1°C at 40°C
Temperature setting	Rotary knob
Temperature scale graduation	2°C
Set point accuracy	±2°C
Max temperature variation between blocks	0.2°C at 40° C
Heat up time 30° C to 37° C	8 minutes
Heat up time 30° C to 56° C	9 minutes
Heat up time 30°C to 100°C	12 minutes
Number of blocks	2


Electrical supply

<i>Voltage</i>	<i>Cycles</i>	<i>Power</i>
220V-240V	50Hz-60Hz	400W
110V-120V	50Hz-60Hz	400W


OPERATION

Preparation

1. Select the insert block(s) appropriate to the application. A list of available blocks is given on page 21 of this manual. Ensure that both the underside of the block(s) and the top of the hotplate are clean; efficient heat conduction between these two surfaces is essential.


Top view of model DB-2P, showing two insert blocks

The DB-1M with a block for Hi-Temp 96 microplates

2. Place the aluminium alloy insert blocks into the well of the unit and place the tubes containing the sample liquid in the blocks.
3. The DB-1M takes a single block for Techne HI-TEMP 96™ microplates. This unit can take other blocks for other makes of plate, but you must ensure that the plates can withstand the maximum temperature of the protocol. (note that the lid is closed during use to improve thermal contact)
4. The heater design, temperature sensor and control circuit give good temperature control and uniformity, but make sure that there is a close fit of the tubes in the block to allow efficient heat transfer.
5. If less than the maximum number of insert blocks are fitted it is recommended that you fit plastic half blocks to the space either side of the aluminium block. A pair of Plastic half blocks are supplied with the DB-2A, DB-2D and the DB-2P. The dummy blocks have the effect of reducing heat losses, thus improving temperature stability.

The DB-2A showing the use of plastic half blocks

6. Plug the mains cable into the socket in the rear of the unit. Connect to the mains electricity supply with the plug provided or one wired correctly for your supply. Switch the power ON from the switch located next to the mains input lead. The mains indicator on the front panel will light (as well as the fuse neon on the mains input unit).

The Overtemperature Cut-out

In the event of a fault occurring which could cause the heater to remain on, power is removed from the heater when the block temperature reaches approximately 110°C.


Setting the Pre-set Temperature

As supplied, all units have presets set to a minimum temperature. In every case, the presets must be set before use.

DB-1M and DB-2P only

To set the presets, follow this procedure:

- Select the preset (1, 2 or 3) to be adjusted with the selector knob
- Place the sample tubes in the block
- Place a thermometer in the sample liquid contained in the tube or directly in the thermometer hole in the block
- Lift and turn the captive plastic cover from the preset to be adjusted, as illustrated below:

Preset Location

- Insert the screwdriver as supplied into the hole revealed by removing the plastic cover
- Turn the screwdriver clockwise to increase the preset temperature, or anti-clockwise to decrease the preset temperature.

A guide to the number of turns required to attain a particular temperature is given by the graph on page 20. Note that, for units up to Serial N°/4, the temperature range corresponds to approximately 25 complete turns of the potentiometer. For units with Serial N°/5 on, the temperature range corresponds to approximately 4 complete turns of the potentiometer.

- If the setting required is lower than the current setting, use the graph to estimate the number of turns required to reach 10°C *below* the target temperature.
- Allow the temperature of the unit to stabilise after each adjustment.


Remember that the unit heats up quicker than it cools down since the unit is well insulated and heat losses are small. Temperature stabilisation can take approximately thirty to forty minutes for each adjustment.

- Use the graph to estimate the number of turns required to increase the temperature half way to the required temperature.

For example:	up to Ser N°/4	Ser N°/5 on
Current temperature	45°C	45°C
Required temperature	68°C	68°C
From graph number of turn	9	1.2
Turn the preset (turns)	4 1/2	0.6

- Measure the temperature when it stabilises and repeat the step above until the required temperature is attained.

When the required temperature is achieved with sufficient accuracy, replace the cover and make a note of the temperature at which the preset is now set.

1. DB-1M and DB-2P

Select the preset temperature required (1, 2 or 3) with the selector knob. Setting the temperatures on these presets is described on pages 19-20 of this manual.

DB-2A

Turn the temperature control knob to the required temperature. Note that the actual temperature of the block critically depends on the type of block used and the quality of the thermal contact between the base of the block and the unit.

2. The heater (and heater indicator) comes on if the set temperature is higher than the current block temperature.
3. When the set temperature is approached, the heater indicator will begin to flash and when the temperature is stabilised the indicator will flash slowly.
4. Due to variations in heat losses with different designs of insert block, the actual temperature may vary. The units are calibrated using model F3506 insert blocks (or an F4467 in the case of the DB-1M). Running a DB-2A, a DB-2D or a DB-2P with only one insert block and without the plastic half blocks will reduce the maximum achievable temperature by a few degrees.

If you need to control the temperature to a greater accuracy than the instrument's set point, place a thermometer in the special hole in one of the blocks. For greater accuracy still, you can place the thermometer in the sample liquid in one of the test tubes. In either of these cases, it may be necessary to readjust the temperature control knob or preset adjustment to achieve the precise temperature required. Allow the temperature of the unit to stabilise after each adjustment.

5. There will be a time lag between the heater platten and the insert block achieving the set temperature due to thermal contact between them. Guide times are as shown in the Specifications on pages 12 to 14.

After Use

1. When you have finished heating samples, remember that parts of the unit – the tubes, blocks and associated accessories – may be very hot. Take the precautions listed earlier. We recommend that the blocks should be allowed to cool to 70°C before being removed from the Dri-Block®. They will still have to be handled with care
2. If you need to remove an insert block while it is hot, screw the extractor tool into the thermometer hole and lift the block out vertically. Never leave the extractor tool in the block while it is being used in the Dri-Block®.

Removing an insert block with the extractor tool

ADDITIONAL INFORMATION**Insert Blocks**

Insert blocks are made of aluminium alloy and must be ordered separately from the heater units. The following blocks are available from Techne:

<i>Part</i>	<i>Tube Size</i>	<i>Number of Holes</i>	<i>Size</i>
F3501	Plain block	None	95 x 76 x 51
F3502	6mm	30	95 x 76 x 51
F3503	10mm	20	95 x 76 x 51
F3504	12mm	20	95 x 76 x 51
F3505	13mm	20	95 x 76 x 51
F3506	15mm	12	95 x 76 x 51
F3507	16mm	12	95 x 76 x 51
F3508	19mm	8	95 x 76 x 51
F3509	25mm	6	95 x 76 x 51
F3510	10mm cuvettes	2 channels	95 x 76 x 51
F4460	Plain block	Thermometer hole only	95 x 76 x 51
F4461	9mm/7mm	10/20	95 x 76 x 51
F4462	24mm	6	95 x 76 x 51
F4463	26mm	6	95 x 76 x 51
F4464	1.5ml Eppendorf	20	95 x 76 x 51
F4465	0.5ml Eppendorf	30	95 x 76 x 51
F4466	Plastic half block	Plain	95 x 37 x 50
F4467#	HI-TEMP 96™	96 well	95 x 151 x 50
F4468#	Falcon 96 round btm.	96 well	95 x 151 x 50
F4469#	Falcon 96 flat btm.	96 well	95 x 151 x 50
F4470	2 ml Eppendorf	20	95 x 76 x 51
F4471	0.2ml	72	95 x 76 x 51
F4473#	Titertgene	96	95 x 150 x 61
F4474#	Porvair	96	95 x 150 x 48

DB-1M only

Operator Maintenance

NOTE THAT THIS EQUIPMENT SHOULD ONLY BE DISMANTLED BY PROPERLY TRAINED PERSONNEL. **REMOVING THE SIDE, FRONT OR REAR PANELS EXPOSES POTENTIALLY LETHAL MAINS VOLTAGES.** THERE ARE NO OPERATOR MAINTAINABLE PARTS WITHIN THE EQUIPMENT.

In the unlikely event that you experience any problems with your Dri-Block[®] which cannot easily be remedied, you should contact your supplier and return the unit if necessary. Please include any details of the fault observed and remember to return the unit in its original packing. Techne accept no responsibility for damage to units which are not properly packed for shipping: if in doubt, contact your supplier.

1. Cleaning

Before cleaning your unit ALWAYS disconnect from the power supply and allow to cool below 50° C. Your Dri-Block[®] can be cleaned by wiping with a damp soapy cloth. Care should be exercised to prevent water from running inside the unit. Do not use abrasive cleaners.

Before using any cleaning or decontamination method except those recommended here, the responsible body or operator should check with Techne that the proposed method will not damage the equipment.

2. Fuses

Your unit is protected by two fuses. These should only be changed by suitably qualified personnel. If the fuses blow persistently, a serious fault is indicated and you may need to return the unit to your supplier for repair.

DB-1M	230V	2 x F2A;	120/100V	2 x F4A
DB-2A	230V	2 x F2A;	120/100V	2 x F4A
DB-2P	230V	2 x F2A	120/100V	2 x F4A

Replacement Parts

The following parts may be obtained from your Techne dealer:

<i>Part Number</i>	<i>Description</i>
6101308	Extraction tool
FCABLEUK	Mains cable with 13amp UK plug (5 amp fuse)
FCABLEEU	Mains cable with R/A Schuko plug
FCABLEUS	Mains cable with 3 pin US plug
FMW11	Techne HI-TEMP 96™ well plate, pack of 25
FMW12	HI-TEMP 96™ well cover, pack of 50
F4466	Plastic half block
6101325	Screw Driver for Presets