

**N E W**

## ENTERPRISE MODEL 350 PORTABLE pH METER



- Easy to read with back lit liquid crystal display
- Rugged and robust case, IP67 rating
- Reading hold function
- Auto buffer recognition - 4.00, 7.00, 9.22 and 10.00
- 500 hours battery life
- Protective boot and carry case for field use
- Full range of application specific electrodes and consumables

## SPECIFICATIONS

### pH

Range: -2 to +16.00pH  
Resolution: 0.01pH  
Accuracy:  $\pm 0.02$ pH

Input Impedance:  $> 10^{12}$  ohms

### Temperature

Manual Temp, Comp. Range: 0 to 100°C / 32 to 212°F  
Resolution: 0.1°C / 1°F

Auto Buffer Recognition: 4.00, 7.00, 9.22, 10.00  
Calibration: 1 or 2 point

Power: 2 AA cells  
Size: 200(l) x 80(w) x 60(d) mm  
Weight: 370g

**JENWAY**