

**N E W**

## ENTERPRISE MODEL 370 PORTABLE pH METER



- Easy to read with back lit liquid crystal display showing temperature compensated pH or mV
- Rugged and robust case, IP67 rating
- Reading hold function
- Auto buffer recognition - 4.00, 7.00, 9.22 and 10.00
- 500 hours battery life
- Protective boot and carry case for field use
- Full range of application specific electrodes and consumables

## SPECIFICATIONS

<b>pH</b>	
Range:	-2 to +16.00pH
Resolution/Accuracy:	0.01pH/±0.02pH
<b>mV</b>	
Range:	-1999 to +1999mV
Resolution/Accuracy:	1mV/±1mV
Input Impedance:	>10 <sup>12</sup> ohms
<b>Temperature</b>	
Range:	-10 to 105°C/14 to 221°F
Resolution:	0.1°C/1°F
Accuracy:	±0.5°C/±1°F
ATC Range:	0 to 100°C / 32 to 212°F
Manual Temp, Comp. Range:	0 to 100°C / 32 to 212°F
Auto Buffer Recognition:	4.00, 7.00, 9.22, 10.00
Calibration:	1 or 2 point pH
Power:	2 AA cells
Size:	200(l) x 80(w) x 60(d) mm
Weight:	370g

**JENWAY**