

N E W

ENTERPRISE MODEL 470 PORTABLE CONDUCTIVITY & TDS METER



- Easy to read with back lit liquid crystal display
- Rugged and robust case, IP67 rating
- Reading hold function
- Direct calibration on standard solutions or by direct cell constant entry
- Auto selection of conductivity or TDS range for optimum resolution
- 500 hours battery life
- Protective boot and carry case for field use

SPECIFICATIONS

Conductivity

Ranges: 0-19.99 μ S, 0-199.9 μ S, 0-1999 μ S, 0-19.99mS, 0-199.9mS, 0-1999mS
Resolution: 0.01 μ S/0.1 μ S/1 μ S/0.01mS/0.1mS/(1mS*)
Accuracy: \pm 0.5% \pm 2 digits

TDS

Range: 0-19.99mg/l, 0-199.9mg/l, 0-1999mg/l, 0-19.99g/l, 0-199.9g/l, 0-1999g/l
Resolution: 0.01mg/l, 0.1mg/l, 1mg/l, 0.01g/l, 0.1g/l, (1g/l*)
Accuracy: \pm 0.5% \pm 2 digits

Temperature

Range: -10 to 105°C/14 to 221°F
Resolution: 0.1°C/1°F
Accuracy: \pm 0.5°C/ \pm 1°F
ATC Range: 0 to 100°C / 32 to 212°F
Manual Temp, Comp.Range: 0 to 100°C / 32 to 212°F
Power: 2 AA cells
Size: 200(l) x 80(w) x 60(d) mm
Weight: 370g

(* only with cell constant >5)

JENWAY