

N E W

ENTERPRISE MODEL 350 PORTABLE pH METER



- Easy to read with back lit liquid crystal display
- Rugged and robust case
- Reading hold function
- Auto buffer recognition - 4.00, 7.00, 9.22 and 10.00
- 500 hours battery life
- Protective boot and carry case for field use
- Full range of application specific electrodes and consumables

SPECIFICATIONS

pH	
Range:	-2 to +16.00pH
Resolution:	0.01pH
Accuracy:	±0.02pH
Input Impedance:	>10 ¹² ohms
Temperature	
Manual Temp, Comp. Range:	0 to 100°C / 32 to 212°F
Resolution:	0.1°C / 1°F
Auto Buffer Recognition:	4.00, 7.00, 9.22, 10.00
Calibration:	1 or 2 point
Power:	2 AA cells
Size:	200(l) x 80(w) x 60(d) mm
Weight:	370g

JENWAY