

Media Release

Fully featured 63-Zero software for 6320D and 6305

Now 6320D and 6305 spectrophotometers come standard with 63-Zero software and serial cable!

- ◆ Auto setup and recognition
- ◆ Manipulate data or method in offline mode
- ◆ Download methods – upload results
- ◆ Complete replication of instrument functions
- ◆ Store, print or transfer results
- ◆ Compatible with 6300

The 63-Zero software uses secure, bi-directional communication that automatically recognizes the model and serial number of the instrument.

The software can be used with or without the instrument connected. In offline mode, you are able to use your personal computer to create or edit methods that can be uploaded to the instrument when the instrument is connected to the computer.

Using the 63-Zero software, you can take measurements remotely, enhancing the features of the instrument by the capabilities of your computer.

All results, including scans and graphic displays, can be printed from the PC giving clear and hard copy presentation. Results can be exported to Excel which enables use with spreadsheets, word processing, PowerPoint and many other software packages for reporting and presentation purposes.

REGISTER YOUR 30-DAY TRIAL SOFTWARE

www.techneusa.com; follow the link for “**Product Warranty Registration**” under “**Support**” section.

CONTACT:

Techne Incorporated

3 Terri Lane, Suite 10 Burlington, NJ 08016,
Phone: 800-225-9243
www.techneusa.com